

Product Datasheet

HDAC1 antibody (FITC) (orb223779)

Description	Mouse monoclonal antibody against HDAC1 conjugated to FITC
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, ICC
Immunogen	Synthetic peptide: EEKPEAKGVKEEVKLA, corresponding to C terminal amino acids 467-482 of Human and mouse HDAC1
Target	HDAC1
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q13547 , O09106
Entrez	433759 , 3065 , 297893
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.