



Product Datasheet

Histone H3 antibody (orb223784)


Description	Mouse monoclonal antibody to Histone H3
Species/Host	Mouse
Reactivity	Human, Mouse, Yeast
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IP, WB
Immunogen	Synthetic peptide corresponding to Human Histone H3 (C terminal) conjugated to keyhole limpet haemocyanin. Database link: P68431
Target	Histone H3
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG3 kappa
Clonality	Monoclonal
Uniprot ID	P84243, P84244
Entrez	15078, 15081, 3021, 3020
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at 1:1000 in WB
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Line graph illustrates about the Ag-Ab r...



Western blot analysis of using Histone H...