

Product Datasheet

PP2A alpha (methyl-L309) antibody (HRP) (orb223795)

Description

Mouse monoclonal antibody against

to HRP

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

HRP

Tested

IP, WB

Applications
Immunogen

A KLH-coupled peptide with the sequence: KLH-C- $\ddot{A}$ YA-RRTPDYFL-OMe. The c-terminal sequence (as well as the carboxy-methylation) is absolutely conserved from human, mouse to yeast

Target

PP2A alpha (methyl L309)

Preservatives

Purified IgG conjugated to Horseradish Peroxidase using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.01% Thiomersal

Concentration

1 mg/ml

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID

P63330, P63331, P67775

Entrez

5515, 24672, 19052

Dilution Range

Optimal antibody dilution should be determined by titration.

Expiration Date

12 months from date of receipt.