


Product Datasheet

Protein A Cys (orb223805)

Description	Recombinant Protein A Cys
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Protein A Cys
Preservatives	Sterile Filtered White Lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 98% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	33.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Protein A Cys