

Product Datasheet

Protein A Long Form (orb223807)

Description	Recombinant Protein A Long Form
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Protein A Long Form
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	46.6 kDa
Expiration Date	6 months from date of receipt.