

Product Datasheet

Cys-Protein A/G (orb223812)

Description	Recombinant Cys-Protein A/G
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Cys-Protein A/G
Preservatives	Sterile Filtered White Lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE and RP-HPLC analyses.
Source	E. coli
MW	47.8 kDa
Expiration Date	6 months from date of receipt.