

Product Datasheet

Protein L-Cys His tag (orb223815)

Description	Recombinant Protein L-Cys His tag
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Protein L-Cys His tag
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE analysis
Source	E. coli
MW	41.6 kDa
Expiration Date	6 months from date of receipt.