

Product Datasheet

Cys-Protein A/G/L (orb223821)

Description	Recombinant Cys-Protein A/G/L
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Cys-Protein A/G/L
Preservatives	Sterile Filtered White Lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and SEC-HPLC analyses.
Source	E. coli
MW	89.3 kDa
Expiration Date	6 months from date of receipt.