

Product Datasheet

Human IGF-BP5 protein (orb223836)

Description	Recombinant Human IGF-BP5 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human IGF-BP5
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2 µm filtered concentrated solution in 10 mM Sodium Citrate, pH 3.0.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE and HPLC analyses
MW	28.6 kDa
Expiration Date	6 months from date of receipt.