


Product Datasheet

Human IGF-BP7 protein (orb223837)

Description	Recombinant Human IGF-BP7 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human IGF-BP7
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.6, 150 mM NaCl
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	26.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IGF-BP7 prote...