


Product Datasheet

Human Proinsulin C-Peptide Analogue protein (orb223851)

Description	Recombinant Human Proinsulin C-Peptide Analogue protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human Proinsulin C-Peptide Analogue
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	3.6 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Proinsulin C-...