


Product Datasheet

Rhesus Macaque IL-1 alpha protein (orb223852)

Description	Recombinant Rhesus Macaque IL-1 alpha protein
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Rhesus Macaque IL-1 alpha
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	18.1 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rhesus Macaque IL-1...