


Product Datasheet

Murine IFN-lambda 3 protein (orb223872)

Description	Recombinant Murine IFN-lambda 3 protein
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Murine IFN-lambda 3
Preservatives	Sterile Filtered White Lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 98% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	19.6 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine IFN-lambda 3...