


Product Datasheet

Rat IL-2 protein (orb223880)

Description	Recombinant Rat IL-2 protein
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Rat IL-2
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	15.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat IL-2 protein