


Product Datasheet

Human VCC-1 protein (orb223900)

Description	Recombinant Human VCC-1 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human VCC-1
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	11.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human VCC-1 protein