


Product Datasheet

Human TAFA-2 protein (orb223906)

Description	Recombinant Human TAFA-2 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human TAFA-2
Preservatives	Sterile Filtered White lyophilized (freezedried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDSPAGE and HPLC analyses.
Source	E. coli
MW	11.2 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human TAFA-2 protein...