


Product Datasheet

Human PTH7-84 15N protein (orb223925)

Description	Recombinant Human PTH7-84 15N protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human PTH7-84 15N
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 98% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	8900 Da
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human PTH7-84 15N p...