

Product Datasheet

Human ME2 His tag protein (orb223929)

Description	Recombinant Human ME2 His tag protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human ME2 His tag
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	64.4 kDa
Expiration Date	6 months from date of receipt.