


Product Datasheet

Human TIGAR-TAT protein (orb223930)

Description	Recombinant Human TIGAR-TAT protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human TIGAR-TAT
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	31.7 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human TIGAR-TAT pro...