

Product Datasheet

Human IHH C28II protein (orb223946)

Description	Recombinant Human IHH C28II protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human IHH C28II
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	19.8 kDa
Expiration Date	6 months from date of receipt.