



Product Datasheet

MKK2 (phospho-T394) antibody (orb223988)


Description	Rabbit polyclonal antibody to MKK2 (Phospho-T394)
Species/Host	Rabbit
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T394 of human MKK2 protein. The exact sequence is proprietary.
Target	MAP2K2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P36507
Entrez	5605
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IP: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of MKK2 (Phospho-T394)



Immunohistochemical analysis of MKK2 (Phospho-T394)



Immunofluorescent analysis of MKK2 (Phospho-T394)