



Product Datasheet

NFkB p65 (AcK310) antibody (orb223996)


Description	Rabbit polyclonal antibody to NFkB p65 (AcK310)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IH, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K310 of human NF-kappaB p65 protein. The exact sequence is proprietary.
Target	RELA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q04207 , Q04206
Entrez	19697 , 5970
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500, ChIP: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of NF-kappaB p65 (...)



Immunohistochemical analysis of NF-kappa...



Immunofluorescent analysis of NF-kappaB ...