



Product Datasheet

RPS6 (phospho-S235) antibody (orb223999)

Description	Rabbit polyclonal antibody to RPS6 (Phospho-S235)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S235 of human RPS6 protein. The exact sequence is proprietary.
Target	RPS6
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P62753, P62754, P62755
Entrez	20104, 6194, 29304, 100911372
Dilution Range	WB: 1:500-1:1000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of RPS6 (Phospho-S235)



Immunofluorescent analysis of RPS6 (Phospho-S235)