



Product Datasheet

S6K1 (phospho-S447) antibody (orb224002)


Description	Rabbit polyclonal antibody to S6K1 (Phospho-S447)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S447 of human S6K1 protein. The exact sequence is proprietary.
Target	RPS6KB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8BSK8, P67999, P23443
Entrez	83840, 6198, 72508
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of S6K1 (Phospho-S...



Immunohistochemical analysis of S6K1 (Ph...



Immunofluorescent analysis of S6K1 (Phos...