



Product Datasheet

TSC2 (phospho-T1462) antibody (orb224038)


Description	Rabbit polyclonal antibody to TSC2 (Phospho-T1462)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T1462 of human TSC2 protein. The exact sequence is proprietary.
Target	TSC2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q61037, P49816, P49815
Entrez	24855, 7249
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of TSC2 (Phospho-T1462) antibody.



Immunohistochemical analysis of TSC2 (Phospho-T1462) antibody.



Immunofluorescent analysis of TSC2 (Phospho-T1462) antibody.