



Product Datasheet

NFIA antibody (orb224174)

Description	Rabbit polyclonal antibody against NFIA
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human NFIA. The exact sequence is proprietary.
Target	NFIA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P09414, P70257, P97863, Q12857, P08651, O00712, P70255, Q02780, Q14938
Entrez	4774, 25492, 18029, 18027, 4782, 18028, 4784, 18032, 4781
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NFIA expression...



Immunohistochemical analysis of NFIA sta...