

Product Datasheet

Histone H3 (phospho-T11) antibody (orb224184)

Description

Rabbit polyclonal antibody to Histone H3

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rabbit, Rat, Zebrafish

Conjugation

Unconjugated

Tested

IF, WB

Applications
Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T11 of human Histone H3 protein. The exact sequence is proprietary.

Target

HIST1H3A

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P68431, P68433, Q6LED0, P84244, P84228, Q71DI3, P84245, P84243

Entrez

8352, 8356, 8354, 8968, 15078, 8357, 319152, 8358, 8355, 291159, 100361558, 8350, 15077, 8353, 8351

Dilution Range

WB: 1:500-1:1000, IF/ICC: 1:100-1:500

Expiration Date


12 months from date of receipt.

kDa

A B C

250

70
51
38
28
19
16

Western blot
analysis of Histone
H3 (Pho...

Immunofluorescent
analysis of Histone
H3...