

Product Datasheet

PKR (phospho-T451) Peptide (orb224189)

Description	The peptide is used to block PKR (phospho-T451) antibody (orb223990) reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	EIF2AK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb223990. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q03963, P19525
Entrez	5610, 19106
Expiration Date	12 months from date of receipt.