

Product Datasheet

RPS6 (phospho-S235) Peptide (orb224198)

Description	The peptide is used to block RPS6 (phospho-S235) antibody (orb223999) reactivity
Reactivity	Human, Mouse, Primate, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	RPS6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb223999. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62754, P62753, P62755
Entrez	100911372, 6194, 29304, 20104
Expiration Date	12 months from date of receipt.