

Product Datasheet

Delta-NaCH Peptide (orb224202)

Description	The peptide is used to block Delta-NaCH antibody (orb224003) reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SCNN1D
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224003.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51172
Entrez	6339
Expiration Date	12 months from date of receipt.