

Product Datasheet

TGF Beta 1 peptide (orb224255)

Description	TGF Beta 1 peptide
Reactivity	Human, Porcine, Primate
Tested Applications	Blocking
Target	GDF15
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224056.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z0J7, Q99988
Entrez	23886, 9518
Expiration Date	12 months from date of receipt.