

Product Datasheet

ARMCX2 Peptide (orb224258)

Description	The peptide is used to block ARMCX2 antibody (orb224059) reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ARMCX2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224059. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6A058, Q7L311
Entrez	67416, 9823
Expiration Date	12 months from date of receipt.