

Product Datasheet

LYPLA1 Peptide (orb224268)

Description	The peptide is used to block LYPLA1 antibody (orb224069) reactivity
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	LYPLA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224069. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70470, O75608, P97823
Entrez	25514, 10434, 18777
Expiration Date	12 months from date of receipt.