

Product Datasheet

ATF7 Peptide (orb224275)

Description	The peptide is used to block ATF7 antibody (orb224076) reactivity
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	ATF7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224076. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17544, Q8R0S1
Entrez	223922, 11016
Expiration Date	12 months from date of receipt.