

Product Datasheet

TRAF3IP1 Peptide (orb224288)

Description	TRAF3IP1 Peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TRAF3IP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224089.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8TDR0, Q5XIN3, Q149C2
Entrez	74019, 26146, 363286
Expiration Date	12 months from date of receipt.