

Product Datasheet

HSPA14 Peptide (orb224297)

Description	The peptide is used to block HSPA14 antibody (orb224098) reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	HSPA14
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224098. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6AYB4, Q0VDF9, Q99M31
Entrez	51182, 50497, 307133
Expiration Date	12 months from date of receipt.