

Product Datasheet

KLF1/5/7 Peptide (orb224355)

Description	The peptide is used to block KLF1/5/7 antibody (orb224156) reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	KLF1; KLF5; KLF7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224156. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P46099 , Q13887 , Q9Z0Z7 , Q99JB0 , Q13351 , O75840
Entrez	10661 , 93691 , 688 , 12224 , 8609
Expiration Date	12 months from date of receipt.