

Product Datasheet

CK1 alpha Peptide (orb224359)

Description	The peptide is used to block CK1 alpha antibody (orb224160) reactivity
Reactivity	Human, Mouse, Rabbit, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	CSNK1A1; CSNK1A1L
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224160. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BK63 , Q8N752 , P97633 , P48729
Entrez	122011 , 113927 , 1452
Expiration Date	12 months from date of receipt.