

Product Datasheet

Ah Receptor (phospho-S36) Peptide (orb224363)

Description	The peptide is used to block Ah Receptor (phospho-S36) antibody (orb224164) reactivity
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	AHR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224164. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q75NT5 , P30561 , P35869 , A9YTQ3 , P41738 , Q3U1U7
Entrez	25690 , 11622 , 11624 , 498999 , 57491 , 196
Expiration Date	12 months from date of receipt.