

Product Datasheet

EPHB1/2 (phospho-Y594/604) Peptide (orb224365)

Description	The peptide is used to block EPHB1/2 (phospho-Y594/604) antibody (orb224166) reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	EPHB1; EPHB2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224166.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29323, P09759, Q8CBF3, P54763, P54762
Entrez	2047, 270190, 1969, 24338
Expiration Date	12 months from date of receipt.