

Product Datasheet

Histone H3 Peptide (orb224369)

Description	The peptide is used to block Histone H3 antibody (orb224170) reactivity
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	H3F3A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224170. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P84243, P84245, P84244
Entrez	15078, 15081, 100361558, 3021, 117056, 3020
Expiration Date	12 months from date of receipt.