

Product Datasheet

SMAD2/3 (phospho-T8) Peptide (orb224371)

Description	The peptide is used to block SMAD2/3 (phospho-T8) antibody (orb224172) reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SMAD2; SMAD3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224172.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O70436, P84022, Q15796, Q62432, P84025, Q8BUN5
Entrez	17126, 4087, 17127, 29357, 4088, 25631
Expiration Date	12 months from date of receipt.