

Product Datasheet

MARK Peptide (orb224372)

Description	The peptide is used to block MARK antibody (orb224173) reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MARK1; MARK2; MARK3; MARK4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224173. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P27448 , O08678 , Q8VHF0 , O08679 , Q9P0L2 , Q7KZI7 , Q8VHJ5 , Q96L34 , Q03141 , Q8CIP4 , Q05512
Entrez	60328 , 4139 , 226778 , 2011 , 170577 , 57787 , 13728 , 4140 , 17169 , 232944 , 117016
Expiration Date	12 months from date of receipt.