

Product Datasheet

STAT5 (phospho-Y694/699) Peptide (orb224380)

Description	The peptide is used to block STAT5 (phospho-Y694/699) antibody (orb224181) reactivity
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	STAT5A; STAT5B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224181.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62771, P42229, P42232, P51692, P52632, P42230
Entrez	25126, 6776, 24918, 20851, 20850, 6777
Expiration Date	12 months from date of receipt.