



Product Datasheet

Claudin 2 antibody (orb224434)


Description	Rabbit polyclonal antibody to Claudin 2
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rabbit, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:200, IF/ICC: 1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis on mouse brain lys...



Immunohistochemical staining of human ki...



Immunofluorescence analysis of A549 cell...