


Product Datasheet

DNA Polymerase beta antibody (orb224452)

Description	Rabbit polyclonal antibody to DNA Polymerase beta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 cell ly...