



## Product Datasheet

### FASL antibody (orb224460)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to FASL  |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human, Mouse                        |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | IHC-P, WB                           |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:200 |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |



Western blot analysis of 293 cell lysate...



Immunohistochemical analysis of using FA...