


Product Datasheet

FGFR1 Oncogene Partner antibody (orb224461)

Description	Rabbit polyclonal antibody to FGFR1 Oncogene Partner
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...