


## Product Datasheet

### NRG1 isoform-10 antibody (orb224497)

|                            |                                               |
|----------------------------|-----------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NRG1 isoform-10 |
| <b>Species/Host</b>        | Rabbit                                        |
| <b>Reactivity</b>          | Human, Mouse, Rat                             |
| <b>Conjugation</b>         | Unconjugated                                  |
| <b>Tested Applications</b> | IHC-P, WB                                     |
| <b>Note</b>                | For research use only                         |
| <b>Clonality</b>           | Polyclonal                                    |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:200           |
| <b>Expiration Date</b>     | 12 months from date of receipt.               |



Western blot analysis of SK-OV3 cell lys...