



Product Datasheet

PCAF antibody (orb224528)


Description	Rabbit polyclonal antibody to PCAF
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Application notes	PCAF Functions as a histone acetyltransferase (HAT) to promote transcriptional activation. Has significant histone acetyltransferase activity with core histones (H3 and H4), and also with nucleosome core particles. Inhibits cell-cycle progression and counteracts the mitogenic activity of the adenoviral oncoprotein E1A. In case of HIV-1 infection, it is recruited by the viral protein Tat.
Clonality	Polyclonal
Dilution Range	WB: 1:500-3000, IHC-P: 1:50-200, IF/ICC 1:100-500, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 cell ly...



Immunofluorescence analysis of HUVEC cel...



Immunohistochemical staining of human ki...